

TUBING & CONNECTORS

SPHYGMOMANOMETER TUBES



Ø 3 x 7 mm FEG, 1 meter
Ref.A10257



Coiled tube 3 meters
Ref.A10907

PNEUMATICS TOURNIQUET TUBES



Silicone tube Ø 4,5 x 9,5 mm, 1 meter
Ref.A10647




Ø 4,5 x 9 mm, 1 meter
Ref.ATV459




Silicone tube Ø 3 x 7 mm, 1 meter
Ref.ATS37R



Silicone tube Ø 3 x 7 mm, 1 meter
Ref.ATS37B




Coiled tube 3 meters
Ref.A20013




Coiled tube 3 meters
Ref.A20014

CONNECTORS FOR CUFFS


clip system connection with O-ring



Female coupler online fluted 4.8 mm
Material : ACETAL
Without sealing valve - **Ref.A11012**
With sealing valve - **Ref.A11016**



Male connector with O-ring
Material : ACETAL/EPDM
Ref.A11008



Male end online fluted 4,8 mm
with sealing valve
Material : ACETAL/EPDM
Ref.A11010



O-ring for acetal male connector
Material : EPDM
Ref.A11322



Extension 3 m + male and femelle connectors
Material : SILICONE/ACETAL
Bleu - **Ref.A20001**
Rouge - **Ref.A20006**




Autoclavable, underwater extension 3 m +
male and femelle connectors
Material : SILICONE/ACETAL
Blue - **Ref.A20003**
Red - **Ref.A20005**




Male connector with O-ring
Material : CHROME BRASS
Ref.A11008M

Standard luer lock connection



Luer-Lock female connection
Material : CHROME BRASS/ACETAL
Ref.A11312



Luer-Lock male connection
Material : CHROME BRASS/ACETAL
Ref.A11311




Male / female Luer-Lock connection
Material : CHROME BRASS/ACETAL
Ref.A11354




Extension 3 m + male and femelle connectors
Material : CHROME BRASS/ACETAL/SILICONE
Ref.A20021

Connection with clip system



Female connector with clip system
Material : CHROME BRASS
Ref.A11305



Male connector with clip system
Material : CHROME BRASS
Ref.A11313




Male / female with clip system
Material : CHROME BRASS
Ref.A11352




CONNECTORS FOR CUFFS AND NIBP


Screw connection



Male connector
Material : ACETAL
Ref.A11306



Female connector
Material : CHROME BRASS
Ref.A11307




Male / female connectors
Material : CHROME BRASS/ACETAL
Ref.A11350


Protective caps for Luer Lock connection



For female connector
Material : PVC
Ref.A45000



Female connector + caps
Material : PVC/NYLON
Ref.A11360




Female connector + caps
Material : PVC/CHROME BRASS
Ref.A11361


Quarter-turn connection



Female connector
Material : ACETAL
Without sealing valve - **Ref.A11301**
With sealing valve - **Ref.A11300**




Male connector quarter-turn connection
Material : ACETAL/EPDM
Ref.A11308




Male / female connectors
Material : ACETAL/EPDM
Ref.A11351

Connection luer lock colin




Female connector
Material : CHROME BRASS
Ref.A11051


Connection with clip system




Female connector
Material : CHROME BRASS
Ref.A11305



Male connector
Material : CHROME BRASS
Ref.A11313



Male / Female connectors
Material : CHROME BRASS
Ref.A11352




Male connector
Material : NYLON
Ref.A11304




Male / Female connectors
Material : CHROME BRASS/NYLON
Ref.A11304


Standard Luer Lock connection




Female connector
Material : CHROME BRASS
Ref.A11312




Male connector
Material : CHROME BRASS
Ref.A11311



Male / Female connectors
Material : CHROME BRASS
Ref.A11354



Female connector
Material : NYLON
Ref.A11303




Male connector
Material : NYLON
Ref.A11302




Male / Female connectors
Material : NYLON
Ref.A11355


Luer connection




Male connector
Material : NYLON
Ref.A11050




Male / Female connectors
Material : NYLON
Ref.A11356




Male / Female connectors
Material : CHROME BRASS/NYLON
Ref.A11357



Male connector
Material : CHROME BRASS
Ref.A11400



Male / Female connectors
Material : LAITON CHROMÉ
Ref.A11358



Male / Female connectors
Material : CHROME BRASS/NYLON
Ref.A11359