




PUSHING TECHNOLOGY TO EXCELLENCE

SIMEON MEDICAL YOUR SOLUTION PROVIDER FOR THE HOSPITAL



Sim.MOVE 800

The next generation
of mobile operating tables



S.I.M.E.O.N. Medical GmbH & Co. KG
In Grubenäcker 18
78532 Tuttlingen
Phone: +49 7461 90068-0
Fax: +49 7461 90068-900
www.simeonmedical.com
info@simeonmedical.com



As a hospital solution provider, SIMEON Medical has over 15 years of experience, 100 employees at 6 sites worldwide, extraordinary innovative capacity and invaluable know-how for medical technology solutions. Thanks to our proximity to our customers and understanding of their needs, we guarantee you excellent products – made in Germany.

Art.No. 500-0018348, CP#003-01384
11/2018


Sim.MOVE 800 – BUSINESS-CLASS WITH HIGH-END CHARACTERISTICS

The Sim.MOVE 800 mobile operating table is fully new developed and made by SIMEON Medical. It is, thanks to its first-class product highlights, the perfect match for your extraordinary requirements in the OR. SIMEON's product philosophy emphasizes intuitive operation, safety, quality and optimal hygiene properties. With a wide range of adjustment options, the Sim.MOVE 800 meets the requirements of all established operating disciplines.



Quality

Maximum load capacity up to 454 kg



Durability, reliability and optimum hygiene properties – SIMEON Quality made in Germany

Easy handling

Intuitive manual control unit, optionally available as wireless version



Easy attachment and detachment of modular components



100% integration possibility: prepared for all operating room integration systems and OR.net ready

Adjustment

Optimum access to the patient and simplified imaging due to longitudinal adjustment of the tabletop up to 400 mm




Enormous flexibility due to large adjustment ranges



Extensive support for common surgical interventions



Ideal working conditions due to the tremendous height adjustment range (600–1100 mm) in a single table version



Functions




Easy maneuverability due to optional fifth wheel

Optional traction drive function with speed adjustment



myLIGHT – Clear and unambiguous coded light signals to enhance user friendliness and safety



AFS – Automatic Floorlock System – Hydraulic floor fixation with levelling of uneven floors for increased stability